





K-Tape[®]... Always a Step Ahead

K-Tape® is an innovative product, recognized worldwide for its optimum quality and performance. The Kumbrink developed original K-Tape® combines the best natural, raw materials with the most advanced manufacturing process available.

- Elastic cotton tape for professional use
- Water resistant
- Sweat-permeable & breathable
- Medication free
- No active ingredients
- Physiobond® acrylic glue
- 96% cotton, 4% spandex





High-End Quality

The Cotton

Only high-quality cotton from one supplier is used for K-Tape. This quality guarantees a fabric structure with homogeneous thread density and virtually knot-free manufacture.

The Elastic Thread

If a thread loses its elasticity, the therapeutic effect is lost. This is why we use the best quality elastane thread for K-Tape. This ensures consistent elasticity levels throughout the entire period of use.

The Dyes

You need to be certain that dyes and adhesive do not harbor risks to the skin or health when applying a tape adhering to the skin for up to 7 days. We use high-quality dyes for *K-Tape*, supplied by Dystar. The reason: Dystar products possess the *OEKO-TEX Standard 100 label* and are safety approved for use in baby clothing. K-Tape is suitable for a variety of applications and a wide range of users.



MDR 2017/745 Ready

The Physiobond® Adhesive

Only *K-Tape* is equipped with professional grade Physiobond adhesive.

Physiobond is a registered trademark in Europe and the USA and is used exclusively for *K-Tape*. We recognize the quality of an adhesive depends not only on the choice of raw materials, but also in the manufacturing and finishing processes.

The objective during production is to create long polymer chains enabling the adhesive to anchor effectively onto the skin. This is followed by a complex temperature-time procedure to remove skin-irritating residual monomers. It is only then, that you can create a first-class tape that establishes a new standard.



Typically, backing paper is only there to protect the adhesive coating. However, with *K-Tape*, the backing paper has additional practical functions as well.

Paper Fibers

The correct orientation of paper fibers makes the backing paper more stable and *K-Tape* pre-cuts deform significantly less than those with conventional backing paper.

Prestretch

It is necessary to prestretch the cotton fabric during production to ensure crease-free manufacturing. It's an important production procedure we use for quality assurance.



K-Tape is applied to the backing paper with 10% prestretch, automatically giving it the ideal amount of prestretch for muscle technique taping. By using the correct technique to apply the tape, it is possible to transfer the required amount of prestretch to the skin. This enables K-Tape to always be pulled off the backing paper with the same amount of tension. Its handling always feels the same... every time you use it.

The Cut

Quite simply, it is the cut that ultimately demonstrates the outstanding quality of *K-Tape*.

A clean, straight cut not only looks professional, it also guarantees better adhesion and duration of wear. If the threads are torn by an unclean cut or if a skewed position causes the outer threads to be frequently cut, the tape unravels at a faster rate and loses its adhesiveness. We therefore take time in production, following a meticulous approach and using high-quality tools, to produce the best cut *K-Tape*



Patients receive a high-quality product, recognized worldwide for its performance quality. Skin irritation is nearly non-existant, which means no patient complaints. Its natural materials contain no harmful substances whatsoever.

K-Tape may potentially adhere to the skin for one week. Many inexpensive tape substitutes adhere for just 2 to 3 days, meaning that they have to be changed frequently. Since patients usually have to pay for all treatments themselves, frequent renewal means higher costs for both patient and therapist.

K-Tape Red

Box of 4 Rolls (20m) - Article # 100141

K-Tape Blue

Single Roll (5m) - Article # 100112



K-Tape Red Single Roll (5m) - Article # 100111



K-Tape Blue Box of 4 Rolls (20m) - Article # 100142





The Original

K-Tape Beige

Single Roll (5m) - Article # 100113

10



K-Tape Black

Single Roll (5m) - Article # 100114



K-Tape Beige

Box of 4 Rolls (20m) - Article # 100143



K-Tape Black

Box of 4 Rolls (20m) - Article # 100144



K-Tape Blue XXL

Single Roll (22m) - Article # 100162



K-Tape Red XXL

Single Roll (22m) - Article # 100161



K-Tape Beige XXL

12

Single Roll (22m) - Article # 100163



K-Tape Black XXL

Single Roll (22m) - Article # 100164





K-Tape Mixed Colors

Box of 4 Rolls (5m each) - Article # 100145







K-Tape My Skin Beige

Single Roll (5m) - Article # 100115



K-Tape My Skin Light Brown

Single Roll (5m) - Article # 100116



K-Tape My Skin Beige

Box of 4 Rolls (20m) - Article # 100146



K-Tape My Skin Light Brown

Box of 4 Rolls (20m) - Article # 100147



K-Tape My Skin Medium Brown

Single Roll (5m) - Article # 100117



Invisible Natural Colors

K-Tape My Skin Medium Brown

Box of 4 Rolls (20m) - Article # 100148



19

K-Tape My Skin Dark Brown

Single Roll (5m) - Article # 100118



K-Tape My Skin Very Dark Brown

Single Roll (5m) - Article # 100119



K-Tape My Skin Dark Brown

Box of 4 Rolls (20m) - Article # 100149



K-Tape My Skin Very Dark Brown

Box of 4 Rolls (20m) - Article # 100150



K-Tape My Skin Beige XXL

Single Roll (22m) - Article # 100176



Single Roll (22m) - Article # 100177





K-Tape My Skin Medium Brown XXL

Single Roll (22m) - Article # 100178



K-Tape My Skin Dark Brown XXL

Single Roll (22m) - Article # 100179

20



K-Tape My Skin Very Dark Brown XXL

Single Roll (22m) - Article # 100180





K-Tape My Skin Mixed Colors

Box of 5 Rolls (5m each) - Article # 100151



21



25

Feel the Benefits of Natural Dye-free Cotton

Since its introduction. K-Tape Pure has gained considerable attention and success, as many of our customers tried the innovative new oxygen-bleached tape right away. Their reaction and approval were immediate and we received extraordinary positive feedback. They discovered its three new distinct benefits, which are a major first in the industry.

Softer Touch

K-Tape Pure feels softer than our well-known original K-Tape. The reason? Its woven pure cotton fabric is not color-dyed. Made from tiny cotton fibers spun into a thread, the resulting tiny spaces its weave creates gives it natural elasticity and a noticeable softer feel. If the thread were dyed, those spaces would fill up with color and the fabric would lose some of its natural flexibility and feel firmer. By being 100% dye-free, K-Tape Pure retains its natural flexibility and is softer to the touch. It's a significant difference that no other tape in the industry can deliver.



K-Tape PureSingle Roll (5m) - Article # 100120

Stronger Adhesion

K-Tape Pure adheres even better, yet still feels light on the skin. It's made with our proprietary Physiobond adhesive, like all our K-Tape products. However, its difference is in its undyed cotton fabric. This results in an improved adhesive quality and establishes an optimal connection between the adhesive and cotton. The open-pore structure of the undyed cotton in K-Tape Pure, provides a key advantage with its unique innovative weave. The adhesive is able to penetrate much deeper into the open spaces and anchor there much better, compared with colour-dyed tapes.

Longer Connection

Patients report high body-tape feedback that describes the lasting effectiveness of K-Tape Pure, even after several days. Typically, the mechanical influences from exercise, normal physical activity and showering impact tape adhesion performance. This can cause the adhesive to partially detach from the dyed-smooth thread cotton surface. However, since adhesive can better penetrate K-Tape Pure's open-space, pure cotton weave, it interlocks more firmly. The better the adhesive anchors in the cotton, the longer the skin connection lasts. This results in a longer-lasting performance with K-Tape Pure.



K-Tape PureBox of 4 Rolls (20m) - Article # 100152




K-Tape Pure XXLSingle Roll (22m) - Article # 100181



Innovative Therapy Treatments

The Kumbrink-created Crosstape patch is used to stimulate selected acupuncture points, enabling an acupuncture-like treatment without the use of needles. It can effectively be used on trigger points, acupuncture points, tense muscles, sore joints and more.







How Does Crosstape® and Crosspuncture® Therapy Work?

The Kumbrink-developed *Crosstape* patch is made of dual mixed fibers and is attached to a specially coated backing paper. When pulled off the backing paper, the patch is charged up electrostatically and ready for specific *Crosspuncture Therapy*. After, it is applied to the skin to stimulate acupuncture points and the areas of pain, where it slowly releases its electrical charge. *Crosstape* patches are waterresistant, can be worn for several days and do not contain any active ingredients.


Many Sizes. Many Applications.


The Crosstape patch is available in sizes S, M, L and XL. Size M is used most of the time.

Size S is ideal for facial applications, acupuncture points on fingertips, tips of toes and ears.

Sizes L and XL can be used effectively on larger areas.

(The 4 Crosstape patches illustrated below are shown actual size.)





Crosstape Mix

Box of 6 Sheets (55 Patches) - Article # 200224

Discover Benefits, Usage and Diverse Applications

For detailed information on how to properly use *Crosstape* patches and how *Crosspuncture Therapy* works to effectively deliver a variety of benefits, visit...

www.crosspuncture.com

Common Crosstape Applications

For headache in the temples

Step 1

Step 2



Headaches arising from tension in the neck



For Stomach Problems



Epicondylitis (Tennis Elbow)

Step 1

Step 2



Professional Accessories. Professional Results.



Using professional tools, yields optimum results. Our exclusive K-Scissors have Teflon-coated blades to prevent the penetration of the acrylic adhesive into the cutting edges and to retain long-lasting sharpness. Our PreKGel® delivers superior skin preparation when applying K-Tape®.



K-Scissors K210Special scissors for right-handers
Article # 300101



K-Scissors K210L Special scissors for left-handers Article # 300103



K-Scissors K160n Bandage scissors Article # 300102

PreKGel Stops the Sweat

PreKGel is a revolutionary product created specifically to stop sweating temporarily. When your patient is in a hot area, or an athlete has finished a warm-up, you are unable to apply tape because of wet, sweaty skin. Simply cleaning with alcohol does not help. PreKGel's special formula delivers major benefits. Its performance-driven blend uses alcohol to clean, green tea and lemon to stop the sweat for a few minutes and aloe-vera to moisturize the skin.



PreKGelGel
Article # 600300


35

34


Best-Selling Guide

Birgit Kumbrink
Author, physiotherapist,
founder and medical/
scientific director of the
international K-Taping

Birgit Kumbrink's best-selling guide to K-Taping® has become an industry standard. These richly-illustrated guide features highly successful therapeutic approaches to treating orthopedic, traumatolgical and many other conditions. It serves as a



Your Go-To Reference Sources



K-Taping in PediatricsBy Birgit Kumbrink (1st edition)
Article # 603104

